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Dynamics of chromosome spreading

&

## Chromprep S4 恒温系统

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Copyright 2012 Sinochrome. All rights reserved

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

### Criteria for Spreading Chromosome

- ✓ 展开面积/Metaphase area
- ✓ 染色体重叠数/Number of Chromosome Overlaps
- ✓ 断裂染色体数/Frequency of Broken Metaphase
- ✓ 带显分辨率/Bands per Haploid Set (BPHS)

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## 如何评价影响因素 How to evaluate these factors

U W Grummt *friedrich schiller university, Germany*

It is well known to experienced cytogeneticists that the quality of metaphase spreads is mainly determined by evaporation of the fixative from the slide, a process that is influenced by climate conditions, such as **relative humidity(湿度) and room temperature(温度)**.

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*Cytogenet Genome Rcs 98:136-146(2002)*

Wen Deng *University of Hong Kong China*

**Humidity (湿度)** over the slide was the most important variable affecting the quality of chromosome spreads.

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*Cytometry Part A 51 A:46-51 (2003)*

Jack L. Spurbeck

*Division of Laboratory Genetics and Section of Biostatistics,  
Mayo Clinic and Mayo Foundation, American*

Optimum metaphase areas can be achieved at various combination of relative **humidity** and **temperature**. We propose that the use of an environmental control unit is a practical way of achieving optimum chromosome spreading routinely and in a highly consistent manner.

*American Journal of Medical Genetics 61:387-393 (1996)*

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
Metaphases are sometimes incomplete with chromosome loss or chromosomes are embedded into the cytoplasmatic matrix with poor banding quality. Such a variability is a common problem in laboratories without appropriate air conditioning during the weather changing.

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

Dry at optimum rate and metaphases have few overlaps  
and not broken


固定液蒸发太慢，染色体断裂、卷曲，分得太散  
Dry too slowly and metaphases are broken and rolled.

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
### Multiple Dynamic Steps in Metaphase spreading

- ✓ 0~10 s, 细胞漂移, 进而固定在玻片上  
Cells were seem to float and then touch the slide (immobile)
- ✓ 10~25 s, 染色体无可见展开  
No visible changes in dimensions of metaphases and chromosomes do not spread.
- ✓ 25~50 s, 染色体展开, 最后5s, 展开最快  
Dimensions of cells increase and chromosomes spread
- ✓ 50~90 s, 染色体变得可见  
Chromosomes are visible on the slides.

*From: U. Claussen, Cytogenetic and Genome Research 2002;98:136-146.*

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
Fig. 7. Estimated median metaphase area as a function of temperature (°C) and relative humidity (%) for lymphocytes based on proportional hazards regression analysis.

*From: American Journal of Medical Genetics 61:387-393 (1996)*

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Fig. 5. Estimated median metaphase area for amniocytes as a function of temperature (°C) and relative humidity (%) based on proportional hazards regression analysis.

*From: American Journal of Medical Genetics 61:387-393 (1996)*

## 染色体断裂—温度+湿度 (羊水细胞) Breakens—Humidity&Temperature(Amniocytes)


Relative Humidity(%)	Temperature		
	20	25	30
25			0/90
26		0/30	
30			
35	0/180	0/150	3/150
45	3/150	13/150	16/150
55	11/150	21/150	22/150
65	26/150	26/150	38/150

*From: American Journal of Medical Genetics 61:387-393 (1996)*

## 染色体断裂—温度+湿度 (淋巴细胞) Breakens—Humidity&Temperature (lymphocytes)

Relative Humidity(%)	Temperature		
	20	25	30
35			2
38		0	
45	0	1	4
55	3	2	13
65	0	8	22


*From: American Journal of Medical Genetics 61:387-393 (1996)*

## Influence of temperature and humidity on the dynamics of chromosome spreading

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


## 国内相关医院状况 Conditions of Chinese Related Hospitals

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## 参考文献 References

1. Jack L., Dynamics of chromosome spreading. American Journal of Medical Genetics 1996;61:387-393.
2. AST Lim, The importance of high resolution chromosome analysis in the diagnosis of birth defects: case reports of holoprosencephaly and cystic hygroma. Annals Academy of Medicine 2004;33:537-540.
3. Octavian Henegariu, Improvements in cytogenetic slide preparation: controlled chromosome spreading, chemical aging and gradual denaturing. Cytometry 2001; 43:101-109.
4. U. Claussen, Demystifying chromosome preparation and the implications for the concept of chromosome condensation during mitosis. Cytogenetic and Genome Research 2002;98:136-146.

- 5 **Wen Deng, A new method for improving metaphase chromosome spreading. Cytometry Part A 2003;51A:46-51.**
- 6 **6R. Hliscs, The spreading of metaphases is a slow process which leads to a stretching of chromosomes. Cytogenetics and Cell Genetics 1997;76:167-171.**
- 7 **A. C. Van Prooijen-Knegt, Spreading and staining of human metaphase chromosome on aminoalkylsilane-treated glass slides. Histochemical Journal 1982;14:333-344**

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