

# Chromed PR 酚红溶液



**Chromed PR 酚红**用作pH指示剂。pH值低于或等于6.4的酚红溶液显示黄色，pH值高于或等于8.2的酚红溶液显示红色。细胞培养液中常常加入酸红，用于指示培养液pH值的变化，常规使用浓度为11 mg/L。

## 『理化性质』

- 分子量: 354.4
- MF: C<sub>19</sub>H<sub>14</sub>O<sub>5</sub>S
- CAS NO.: 143-74-8

## 『溶液组成』

编号	组分	分量g/L
1	Phenol Red	5
2	KCl	0.2
3	KH <sub>2</sub> PO <sub>4</sub>	0.2
4	NaCl	8
5	Na <sub>2</sub> HPO <sub>4</sub>	2.16

## 『保存』

- 常温保存

## 『规格』

- 无菌过滤
- 浓度0.5%
- 包装规格: 100mL
- 货号CM-05-100

## 『注意事项』

- 仅供科研用途，不作为药物、家庭备用或其它用途。

## 『参考文献』

1. J. Biol. Chem., 236, 589 (1961).
2. The Sigma-Aldrich Handbook of Stains, Dyes & Indicators, Green, F. J., ed., Aldrich Chemical Co.(Milwaukee, WI: 1990), p. 565-567.
3. The Merck Index, 12th ed., Entry# 7397.
4. Berthois, Y., et al., Phenol red in tissue culture media is a weak estrogen: implications concerning the study of estrogen-responsive cells in culture. Proc. Natl. Acad. Sci. USA, 83(8), 2496-2500(1986).
5. Bindal, R. D., et al., Lipophilic impurities, not phenolsulfonphthalein, account for the estrogenic activity in commercial preparations of phenol red. J. Steroid. Biochem., 31(3), 287-293 (1988).
6. Pick, E., and Mizel, D., Rapid microassays for the measurement of superoxide and hydrogen peroxide production by macrophages in culture using an automatic enzyme immunoassay reader. J. Immunol. Methods, 46(2), 211-226 (1981).